Clean-Trace® Test

-  (EN) Surface Protein Plus Test
-  (FR) Test Surface Protein Plus
-  (DE) Oberflächen Protein Plus Test
-  (IT) Test Surface Protein Plus
-  (ES) Test de Proteína para Superficie Plus
-  (NL) Oppervlakte eiwittest plus
-  (SV) Surface Protein Plus Test
-  (DA) Overflade Protein Plus Test
-  (NO) Proteintest til overflate
-  (FI) Surface Protein Plus Pintatesti
-  (PT) Teste de Proteína Plus de Superfície
-  (EL) Τεστ Επιφανειών Protein Plus
-  (PL) Test do wykrywania białek na powierzchni „Plus”
-  (HU) Felszíni Fehérje Plusz Teszt
-  (CS) Test Povrchový protein Plus
-  (RO) Proteine de suprafață Test Plus
-  (RU) Тест для определения белка на поверхности (протеин плюс)
-  (TR) Plus Yüzey Protein Testi
-  (JA) タンパク残留測定
-  (ZH) 表面蛋白质 Plus 采样棒

Clean-Trace® Test

Product Instructions

Surface Protein Plus Test

Product Description and Intended Use

Neogen® Clean-Trace® Surface Protein High Sensitivity Test detects protein residues. The results of the tests are time and temperature dependent. By adjusting the time and / or temperature it is possible to achieve different limits of detection. Selecting the correct incubation time and temperature can therefore optimize the sensitivity of the test. A 45 minute incubation at 37°C (99°F) results in equivalent sensitivity to the Ninhydrin test specified in the UK HTM2030 guidelines. If your requirements for sensitivity are different please use the Time Temperature matrix below as a guide.

Temperature / Time (Minutes) / Sensitivity (micrograms)				
	10 Minutes	15 Minutes	30 Minutes	45 Minutes
20°C (68°F)	50 µg	30 µg	20 µg	10 µg
37°C (99°F)	20 µg	10 µg	5 µg	3 µg

Based on the Biuret reaction. Under alkaline conditions the copper ions (Cu^{2+}) form a complex with peptide bonds of protein and become reduced to Cu^+ . Bicinchoninic acid (BCA) under alkaline conditions is a highly sensitive, stable and specific reagent for Cu^+ forming a purple complex. The test provides an estimation of surface protein residue by means of a color change result.

The color may form on the swab more than the solution with certain types of soil. Record the strongest color change on the swab or in the solution as the end result.

Neogen Clean-Trace Surface Protein High Sensitivity Test can be used to assess the level of contaminating protein residue on instruments after they have been processed in a washer/ disinfectant. The test should be used after cleaning to assess whether instruments are clean enough for effective sterilization.

⚠ Warnings

Consult the Safety Data Sheet for additional information and local regulations for disposal.

⚠ Cautions

Protein is a common substance. It is present on skin, hair and on many surfaces. Do not touch the swab or the surface to be tested.






Do not use Neogen Clean-Trace Surface Protein High Sensitivity Test past expiry date.

Although the Neogen Clean-Trace Surface Protein High Sensitivity Test is very sensitive, a surface cannot be considered sterile on the basis of a Neogen Clean-Trace Surface Protein-High Sensitivity Test result.

There is minimal interference with most common sanitizers at normal concentrations. However high concentrations of alkaline cleaners may cause false negative results in some instances. Peroxide based disinfectants can cause a false positive color change reaction with the Neogen Clean-Trace Surface Protein High Sensitivity Test.

User Responsibility

Users are responsible for familiarizing themselves with product instructions and information. Visit our website at www.neogen.com, or contact your local Neogen representative or distributor for more information.

When selecting a test method, it is important to recognize that external factors such as sampling methods, testing protocols, sample preparation, handling, and laboratory technique may influence results.

It is the user's responsibility in selecting any test method or product to evaluate a sufficient number of samples with the appropriate matrices and microbial challenges to satisfy the user that the chosen test method meets the user's criteria.

It is also the user's responsibility to determine that any test methods and results meet its customers' and suppliers' requirements.

As with any test method, results obtained from use of any Neogen Food Safety product do not constitute a guarantee of the quality of the matrices or processes tested.

Limitation of Warranties / Limited Remedy

EXCEPT AS EXPRESSLY STATED IN A LIMITED WARRANTY SECTION OF INDIVIDUAL PRODUCT PACKAGING, NEOGEN DISCLAIMS ALL EXPRESS AND IMPLIED WARRANTIES, INCLUDING BUT NOT LIMITED TO, ANY WARRANTIES OF MERCHANTABILITY OR FITNESS FOR A PARTICULAR USE. If any Neogen Food Safety Product is defective, Neogen or its authorized distributor will, at its option, replace or refund the purchase price of the product. These are your exclusive remedies. You must promptly notify Neogen within sixty days of discovery of any suspected defects in a product and return it to Neogen. Please contact your Neogen representative or authorized Neogen distributor for any further questions.

Limitation of Neogen Liability

NEOGEN WILL NOT BE LIABLE FOR ANY LOSS OR DAMAGES, WHETHER DIRECT, INDIRECT, SPECIAL, INCIDENTAL OR CONSEQUENTIAL DAMAGES, INCLUDING BUT NOT LIMITED TO LOST PROFITS. In no event shall Neogen's liability under any legal theory exceed the purchase price of the product alleged to be defective.

Storage and Disposal

For maximum shelf life, store between 2-25°C (36-77°F). Alternatively, swabs are stable for 14 days when shipped or stored at elevated temperature; for this storage option do not exceed 35°C (95°F). Expiry date and lot number are noted on each pouch. After use, may contain microorganisms that may be a potential biohazard. Follow current industry standards for disposal.

Instructions for Use

Before beginning the test procedure, ensure the dry bath, water bath, or incubator is turned on and set for 37°C (99°F).

1. Place the Neogen Clean-Trace Surface Protein High Sensitivity Test at room temperature for at least 10 minutes before use. Grip the sample stick handle and remove the swab from the device. Use the moisturizer provided and apply 4 drops to the swab OR apply 4 drops to the surface to be tested [figure 1].
2. Swab the test instrument. Where practical, swab an area approximately 10 cm by 10 cm and swab the area in one direction and then swab the area in the opposite direction [figure 2]. Test hard to reach places like serrated edges and joints. Apply pressure on the swab and rotate the swab as the sample is collected to ensure repeatable and effective sampling [figure 3].
3. Reinsert the sample stick into the test with the handle inserted to the original position of the unused test. At this point you can either activate and measure immediately, or, if it is more convenient, the Neogen Clean-Trace Surface Protein High Sensitivity Test may be labeled and left un-activated for up to four hours before activation and measurement of multiple swabs simultaneously.
4. To process the sample, push down firmly on the top of the Neogen Clean-Trace Surface Protein High Sensitivity Test sample stick handle [figure 4]. The handle will slide into the device tube and the top of the handle should be level with the top of the device tube when fully depressed. Grip the top of the device and shake rapidly side-to-side for at least five seconds to mix the sample and reagent [figure 5]. The solution will turn mint green.
5. Place the Neogen Clean-Trace Surface Protein High Sensitivity Test into a heating block that has been pre-heated to 37°C (99°F). Heat for 45 minutes at 37°C (99°F). Compare the color of the swab and the solution against the Neogen Clean-Trace Surface Protein High Sensitivity Test label. Record the strongest color change in the solution or swab as the end result.

How to best utilize Neogen Clean-Trace Surface Protein High Sensitivity Test	
Do	Do Not
Hold the Neogen Clean-Trace Surface Protein High Sensitivity Test in the VERTICAL position on activation and shake rapidly from SIDE TO SIDE for at least 5 seconds.	Activate when not in the vertical position, do not shake up and down, do not shake aggressively.
Incubate Neogen Clean-Trace Surface Protein High Sensitivity Test vertically.	Incubate Neogen Clean-Trace Surface Protein High Sensitivity Test horizontally.
Read Neogen Clean-Trace Surface Protein High Sensitivity Test at 45 minutes.	Record any color change after 45 minutes.

Interpretation

The color of the test reaction will indicate the level of protein residues on the surface. By comparing the color produced against the Neogen Clean-Trace Surface Protein High Sensitivity Test label an estimation of the surface cleanliness can be made.

- Green indicates a Pass result - Clean, No further action required
- Grey or Purple indicates a Fail result, Re-clean & Re-test

For further information, contact Neogen Food Safety Technical Services.

Explanation of Symbols

info.neogen.com/symbols

Neogen Food Safety

Neogen Corporation

620 Leshar Place
Lansing, MI 48912 USA
www.neogen.com

Neogen Europe Ltd.

The Dairy School
Auchincruive
Ayr, KA6 5HU
Scotland, UK

Neogen Ireland, Ltd.

Bray Business Park, Bray
Co. Wicklow
A98YV29, Ireland

Neogen Corporation

620 Leshar Place Lansing, MI 48912 USA
www.neogen.com

© Neogen Corporation 2024. All rights reserved.
Neogen and Clean-Trace are registered trademarks of Neogen Corporation.
FS00926A

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